


Uploading CT Scan Data to 3D Gens Server Guideline


Please go to website as following url address:

nas3dgers.quickconnect.to




Please enter the username and password to login :


Username: To be confirmed
Password: To be confirmed




Click at the dedicated file (provided by 3d gens) and click at "CT scan data" file :




Crete new folder via "Create Folder" and rename the file according to the patient name. Enter the patient name file by clicking the file twice:




Click twice at the "File Station" :




Click upload button and choose upload skip:




Choose and drag the desired file to upload to the server.
Please ensure the file is in DICOM format:



To start upload the file, please drag the file into "patient 1" or the newly created patient's name file as showed in the diagram



Click the skip button to proceed the data uploading process:



Data successfully uploaded as showed in the diagram:

